

**BVI**

# Providing Relief for Patients with Dry Eye or Lacrimal Disorders

Parasol<sup>®</sup> | Micro Flow<sup>™</sup>



**Permanent  
Parasol<sup>®</sup>**



**Partial  
Micro Flow<sup>™</sup>**

# Permanent Occlusion

Permanent occlusion is indicated for chronic dry eye conditions. All occluder products are made of medical grade silicone and are reversible.

## Insertion

The hollow nose collapses on insertion, eliminating the need for punctal dilation in most cases. It's straight-forward for you to insert and more comfortable for your patients.

## Straight Forward Sizing

Two sizes, small and medium, fit virtually all patients. X-Small and large sizes are available for customization.

## Retention

Parasol is designed for excellent retention, as demonstrated by an independent clinical study.<sup>1</sup>

## Parasol® Size Comparison Chart



Traditional Plugs

0.2 mm	0.3 mm	0.4 mm	0.5 mm	0.6 mm	0.7 mm	0.8 mm	0.9 mm +
--------	--------	--------	--------	--------	--------	--------	----------

Parasol Plugs

X-SMALL	SMALL	MEDIUM	LARGE
---------	-------	--------	-------




## Parasol® Pre-Loaded Punctal Occluders Sterile

Occluders are packaged in a convenient 1 pair tray, are sterile and pre-loaded on a single-use inserter.



A14-201




Parasol® Pre-Loaded Punctal Occluders Sterile	SKU #
X-Small, 0.25 - 0.35 mm	A14-200
Small, 0.35 - 0.65 mm	A14-201
Medium, 0.6 - 0.85 mm	A14-202
Large, 0.9 mm+	A14-203

<sup>1</sup> C. McCabe, MD, PhD, FACS Plugs Reduce Dry-Eye Symptoms, Improve Vision. November 2009, Review of Ophthalmology.

## Parasol® Quickload Punctal Occluders Sterile

Occluders are sterile and pre-mounted on a stand for Quickload convenience.

QTY: 16/Bx



Parasol® Quickload Punctal Occluders Sterile	SKU #
Small, 0.35 - 0.65 mm	A14-201QL
Medium, 0.6 - 0.85 mm	A14-202QL
Large, 0.9 mm+	A14-203QL
Combo small and medium	A14-205QL

## Parasol® Bulk Punctal Occluders Non-Sterile

Occluders are non-sterile and packaged 2 occluders per poly bag.

QTY: 20/Bx



Parasol® Bulk Punctal Occluders Non-Sterile	SKU #
Small, 0.35 - 0.65 mm	A14-201BN
Medium, 0.6 - 0.85 mm	A14-202BN
Combo small and medium	A14-209BN

# Partial Occlusion

## Micro Flow™ Punctal Occluders

Designed to provide a straight-forward, reversible means of partial occlusion. Micro Flow™ Punctal Occluders support lacrimal drainage care. Suggest to use when total occlusion is not ideal. Occluders are packaged in a convenient 1 pair tray, are sterile and pre-loaded on a single-use inserter.



A14-401



QTY: 2/Bx

Micro Flow™ Pre-Loaded Punctal Occluders Sterile	SKU #
Sterile, Small, 0.4 - 0.55 mm	A14-401
Sterile, Medium, 0.55 - 0.7 mm	A14-402
Sterile, Large, 0.7 - 0.85 mm	A14-403

QTY: 10/Bx

Micro Flow™ Punctal Occluders Non-sterile	SKU #
Non-Sterile, Bulk, Small, 0.4 - 0.55 mm	A20-401BN
Non-Sterile, Bulk, Medium, 0.6 - 0.85 mm	A20-402BN
Non-Sterile, Bulk, Large, 0.9 mm+	A20-403BN



# Single-Use Punctal Dilators and Probes

## Visitec® Punctal Dilator and Probe, Khan

Dual function instrument allows for single handed transition from dilating to probing without release of lateral lid traction.



**QTY:** 3/Bx  
**Brand:** Visitec®

Description	SKU #
25 mm/1 in medium taper, Blunt, Wilder style dilator mated with a size 0 Bowman lacrimal probe, 0.7 mm x 45 mm/22G x 1 3/4 in	585124

## Punctum Dilator, Fine/Medium

0.4 mm x 21.5 mm Tip,  
0.6 mm x 18 mm Tip,  
77 mm Handle




**Qty:** 20/Bx  
**Brand:** Malosa®

Description	SKU #
Punctum Dilator, Fine/Medium	MMSU1181S

## Punctum Dilator/Seeker

0.45 mm x 21.5 mm Tip,  
5 mm Pointed tip,  
18 mm Straight shaft, 77 mm Handle



**Qty:** 20/Bx  
**Brand:** Malosa®

Description	SKU #
Punctum Dilator/Seeker	MMSU1193S



# Other Related Occlusion Options

## Curved Tying Forceps

0.25 mm Tips, 5 mm Curved tying platforms,  
80 mm Handle



Qty: 20/Bx  
Brand: Malosa®

Description	SKU #	SKU # Clips
Curved Tying Forceps	MMSU1210S	MMSU1210CS


## Lacrimal Cannula



Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 1/2 in/0.5 x 13 mm	585068

## Lacrimal Cannula



Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 26G x 1 1/8 in/0.45 x 28 mm	581276



FOR COMPLETE PRODUCT INSTRUCTIONS,  
PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Questions: [bvimedical.com/customer-support](https://bvimedical.com/customer-support)



**BVI**  
bvimedical.com